

How to Create Your Voice Clone

Follow these steps to create your ElevenLabs Professional Voice Clone.

ElevenLabs Video Tutorial :

A detailed video tutorial is available for a visual walk-through on how to create your voice clone. Access the video [here](#).


ElevenLabs Audio Editing Support :

Request audio editing support for your Professional Voice Clone. Access form [here](#).

Part I: Creating a New Voice Clone


After you receive your acceptance email from ElevenLabs, your account will be automatically upgraded to an ElevenLabs Pro account the next time you log in.

- Go to the Website: Go to elevenlabs.io and click **Sign in**.
- Navigate to the creative platform and click the plus sign next to **Voices**.




How to Create Your Voice Clone

Select the Professional Voice Clone tab.




Click on Create New Clone.



Part 2: Uploading or Recording Voice Samples


Add a name, specify the language, and accent.



How to Create Your Voice Clone

Upload Samples


- Select your files (supported formats include MP3, WAV, M4A, etc., and WhatsApp voice recordings).
- Upload up to 100 files; each must be longer than 30 seconds and good quality.



Crucial Requirement: Ensure you have uploaded at least 10 minutes of total recordings before proceeding.

Record Yourself

- Review guidance and select a script (conversational style).
- Click on Start Recording. Each recording must be at least 30 seconds long.




How to Create Your Voice Clone

Part 3: Voice Verification

Automatic Verification (If your voice is unchanged)

- Repeat the random phrase (you have five attempts).



If successful, your professional voice clone will begin creation (ready in 2 to 6 hours).



Manual Verification (If your voice has changed or verification fails)

- Do not press 'start recording.' Click on Help Center.
- Complete the support request form with your email, subject (e.g., unable to complete voice verification), and description of circumstances. The support form is accessible by clicking [here](#).

How to Create Your Voice Clone

- In the Type of Problem dropdown, select Product, then Professional Voice Clone. In Type of Request, select Question.
 - Please include the following text in the description box:
 - “Due to my multiple system atrophy (MSA) diagnosis, I am currently unable to complete the voice verification process as my current voice quality no longer matches the quality of the original recording. Therefore, I request manual verification.”
 - Click Submit. (ElevenLabs support team may take up to 48 hours to respond).
 - Provide Documents: Complete the form emailed by ElevenLabs, detailing the differences between your current voice and the voice you wish to clone. Save this form along with your photo ID and medical documentation of your condition.
- Log back in: [Voices](#) > [My Voices](#). Click the three dots next to the voice clone name, select Finish This Voice, then Next.
- Click the Manual Verification tab. Click Upload files and select all three documents simultaneously.
- Add a supporting statement. Click Request Verification.
- If verified, the voice clone will start building (ready in 2 to 6 hours).



Your Professional Voice Clone is being prepared

This can take **2 to 6 hours** to complete, and you can track the progress in this [voice's details panel](#). In rare cases the voice-tuning process may fail but it will automatically restart. If it fails and doesn't restart, don't hesitate to [contact support](#).

While the voice is being created, you won't be able to use it for the models it's being fine-tuned on.

[Finish and return to Dashboard](#)